

Operating Instructions for

Helios NanoLab 600

Dualbeam Focused Ion Beam & Scanning Electron Microscope

RULES OF USE

- 1) DO NOT TOUCH A CONTROL IF YOU DON'T KNOW EXACTLY WHAT IT DOES.**
- 2) NEVER, EVER FORCE ANYTHING BEYOND FINGER STRENGTH.**
- 3) IF IN DOUBT - ASK FOR HELP.**

Please note: This instruction can only be used by trained users as a reference. It is not intended to be the textbook on DIY purposes for the unauthorized personnel to use the machine without training and direct supervision.

ANFF@ANU

Front note

- The right hand side computer is used to run operating software of Helios, which should be setup and ready to use
- The left hand side computer has the software for Genesis EDS system and TSL EBSD data collection and analysis system
- The computer in the middle is the node connected to the ANFF server. All the data collected on Helios are saved in local director "SharedData/users" at first, and then uploaded to the server through the node by the user
- Usage will be monitored with both login information and system log files

Task 1: Start a FIB session

- Check list:
 - The sample size is suitable for the machine
 - The sample is flat and locked onto a pin stub
 - The sample is dry and clean
 - The surface is conductive and well earthed
 - The silver paint is fully solidified, if applied
 - The sitemap is available for multiple samples mounted on the same stub


- Double click system control software



- Click on Start UI











- Login your user account






- If the user interface is hidden off the scene on the right hand side computer, click on Show UI on the screen,




Task 2: Load sample and get E beam image






- Vent the chamber
 - Go to page 1, 
 - Click on  and confirm the action in pop-up window
 - Wait for ~ 4 min, then try to open the chamber by pulling the handle bar gently
 - Lock the stub onto the stage with the Allen key, and check the overall height with the gauge, make sure the height lower than the **max line**.
 - Pump the chamber
 - Gently close the door
 - Click on 
 - Wait for ~ 5 min, until the vacuum icon  turn green 
 - Select both electron and ion beam images in serial, and click on  for both columns. If the ion source is not “green”, click on .
- 

The image shows a control panel interface with several sections: 'Vacuum' (Pump, Vent), 'System' (Wake Up, Sleep), and 'Column' (Beam On, Source). A red box highlights the 'Beam On' button in the Column section, with a red arrow labeled '3' pointing to it. A red arrow labeled '2' points to the 'Wake Up' button in the System section. A red arrow labeled '1' points to the 'Source' button in the Column section.
- Locate region of interest (ROI)
 - Select E beam image, e.g. top left quad window
 - Click on pause button  to start E beam scan
 - Reduce the magnification to <100x, and have an overview of whole sample
 - Align the working area if necessary
 - Locate a small feature, e.g. dust, near ROI as the reference
 - Achieve a focused image
 - Link Z to WD
 - Go to page 2, 
 - Click on the link icon  on the top of the window

Task 3: Setup Eucentric height and beam coincidence


- ❖ Standard working conditions are all based on the Eucentric height (~4.0 mm)
- ❖ Pt deposition can only be done at the Eucentric height




- Move stage upwards
 - Increase magnification of E beam image to above 1000x
 - Achieve a fine focus and update Z measurement by click on the icon 
 - Input Z of 8 mm, check focus and update Z measurement with the button 
 - Input Z of 6 mm, check focus and update Z with 
 - Input Z of 4.0 mm, check focus and update Z with 
- Title stage
 - Select "Zero beam shift" from the drop down menu of "stage" for both E beam and I beam image respectively
 - Set a small feature in the centre cross of the E beam image.
 - Input title, T of 15 degrees.
 - Make the feature back to the centre with Z adjustment by the mouse.
 - Input T of 30 degrees, and re-adjust Z to centre the feature.
 - Input T of 52 degrees, and re-adjust Z, if necessary.
- Setup beam coincidence
 - Centre a reference feature on E beam image at 3500x
 - Select I beam image in another window, e.g. top right quad window
 - Check if the I beam current is small, then get a snapshot with the icon 
 - Centre the feature with beam shift knobs if necessary, and update the image with snapshot


Task 4: Cover ROI with Pt deposition

- Select an appropriate ion beam current
 - 2-6 pA/um² is recommended for general application
 - Update the image after 10 seconds if changing the aperture – sleep time to stabilise the beam
- Insert GIS needle
 - ❖ **Make sure that eucentric height has been set up correctly before doing this!**
 - ❖ **Make sure there is nothing higher than ROI nearby!!**

- Go to page 3, 
- Locate the “Pt dep” line, double click on the word “cold” in “Heat” column
- Wait until the yellow timing bar turn to the word “Warm”




- Insert the needle by tick the box “in”
- Confirm the action within the pop up window by click “Yes”



- Define a pattern

- Select rectangular pattern
- Draw a box on the I beam image
 - Green color box for deposition
 - Yellow color box for milling
- Change the Application to “Pt dep”



- Start/stop deposition

- Click on the triangle button  to start deposition
- The same set of buttons will allow the user to pause and stop



- When finish, remove the tick within the box under “in” to retreat the GIS needle

Task 5: Mill cross section and record image

- Select an appropriate ion beam current
 - 6.5 - 21nA are recommended for rough milling
 - For the same pattern, smaller beam current leads to sharper images but in longer milling time



- Define a pattern

- Select/activate “Regular cross section” pattern for rough milling
- Draw a box on the ion beam image
- Confirm Application file as “Si”
 - Other material files are also available in the click/drop-down menu

Name	Value
Application	Si
X size	PMMA
Y size	Pt dep
Z size	Pt dep ebeam
DwellTime	Pt dep ebeam
ScanDirection	SCE
Rel. Int. Diam.(%)	Si
Beam	Si2
TotalTime	Si3N4
	Si_small
	W dep
Progress	W dep ebeam

- Update parameters and location of the pattern if necessary

- Start/stop milling process

- Click on triangle button  to start milling
- The same set of buttons will allow the user to pause and stop 

- Clean surface

- Select a smaller ion beam current, 9.7 – 93 pA are recommended
- Update the image
- Select/activate “Cleaning cross section” pattern
- Draw a narrow box on the ion beam image
- Update parameters and location of the pattern if necessary







- Click on triangle button  to start milling
- The same set of buttons will allow the user to pause and stop 

- Monitor the progress with snapshot 

- Take a photo


- Select/activate the quad desired – E beam imaging is in common, but I beam imaging may show extra information for certain materials, such as channelling contrast for polycrystalline metals
- Check magnification, location and fine focus if necessary
- Press F2 -- pre-set scan speed for photo can be modified within Preferences/Scanning

Task 6: End a session


- Return to Page 1, 
- Turn beam off by clicking on “beam on” for both E beam and Ion beam
- Vent the chamber
 - Click on  and confirm the action within the pop-up window
 - Wait for ~4min, then try to open the chamber by pulling the handle bar gently
- Pump the empty chamber
 - Gently close the door
 - Click on 
 - Wait for ~4min, the vacuum icon  should turn green 
- If you are the last user on the day, turn off the ion source
 - Select/activate ion beam image
 - Click on  to turn off the ion source
- Logout your account

Advanced task: TEM lamella sample preparation

- AutoTEM G2 Software
- Preconditions:
 - Stage Tilt: 52 °
 - High voltage: 30 kV
 - Working distance: Eucentric (defined with E Column)
 - Magnification: 1,000x
- Define the membrane
 - Start the AutoTEM G2 Software
 - Press 'Create' to start the recipe definition process.





- Save the project
 - Select the following recipe '12x5 ex-situ lift-out'
 - Verify the current parameters:
 - Deselect 'skip fiducial milling'
 - Select 'Mill Fiducials' to create fiducials for image recognition.
 - Select 'Refine' to enable accurate lamella positioning
 - Select 'Finish' closing the recipe definition process.
- Run the sites
 - Press the 'Run' button in the position list window.




Advanced task: 3D imaging – Auto Slice & View

- Setup Eucentric height
- Locate ROI and cover it with Pt deposition (with the ion beam current of 2-6pA/ μm^2) – defined a thickness less than 0.5 μm (recommended)
- Make a regular cross section with high beam current, e.g. 9.3 - 21nA, at the leading edge of the ROI – a bit deeper and wider than ROI
- Make trenches on both sides of ROI to avoid shading on images
- Clean up the side edges as the reference for 3D reconstruction
- Select an appropriate beam current for ion beam image, e.g. 93pA or higher
- Using cleaning cross section to clean the viewing area
- Tick the tilt correction in automatic mode
- Hide UI and create a folder in SharedData to store the image files
- Start AutoSlice&View G2 (click on the icon on the desktop)
- Direct image save location to the predefined folder, and change the prefix for data set
- Setup E beam image scan parameter, e.g. 9.05s



- Select Focus Adjust from Setup menu




- Beam Shift for Image Centre is optional
- Select E beam image within the same quad of I beam image by clicking




- Optional step - If high resolution imaging is desired, switch to mode 2 with drop down menu

➤ a range of working conditions including limited kV, nA and WD required to activate mode 2




- From the menu bar, select stage/Zero beam shift
- Double click on the image to centre the viewing area for E beam image
- Return to I beam image, and use beam shift knobs to centre the ROI
- Input parameters for the slice (without Enter)
- Click on Show to see the overlay of the pattern on the ion beam image
- Modify the location with beam shift – do not move the patterns
- Capture/print screen or record important parameters for 3D reconstruction
- HFW for the final E beam image
- Image resolution for E beam image
- Slice length y
- Slice number



- Click on Run and work out the total time
- If there is a pop-up warning for over limitation of image shift, reduce the magnification of E beam image and retry
- If everything is OK, go home!

- *Relationship of the Two Columns*



**PLEASE NOT
REMOVE THIS
INSTRUCTION FROM
THE MACHINE ANYTIME!**